



# MRI JET VENTILATOR

- △ **NEW “Standard of Care” for MRI Suite Airway Emergencies**
- △ **“Compromised” Airways at the MRI can now be Jet Ventilated**
- △ **No need to move non-breathing patient from the MRI suite**
- △ **Tested MR Conditional to 3.0 Tesla**
- △ **Fully Customizable to your precise needs regarding both Input Fittings and Hose Lengths**
- △ **Standard p/n 00-325-MRI Jet Ventilator sold complete with:**
  - **O<sub>2</sub> Input DISS Fitting**
  - **High Pressure Hose 5ft Length**
  - **O<sub>2</sub> Throttling Volume Control Valve**
  - **O<sub>2</sub> Thumb Control Valve ON-OFF**
  - **Back Pressure Gauge 0-60 psi**
  - **Low Pressure Hose 2ft Length**
  - **O<sub>2</sub> Delivery Malleable Copper Stylet 20cm Length**



Brochure 00-325-MRI.doc

## **Anesthesia Associates, Inc.**

460 Enterprise Street  
San Marcos, CA 92078  
Tel: (760) 744-6561 Fax: (760) 744-0054  
Email: Solutions@AincA.com Web: www.AincA.com